



## **BA 904 PORTABLE SLIT LAMP**

High-quality examinations, anywhere...

02 | 03 POWERFUL, HAND-HELD SLIT LAMP

## BA 904

# Haag-Streit optics included

The BA 904 is the only hand-held slit lamp that boasts Haag-Streit's unsurpassed optical system.

The robust, ergonomically-designed portable BA 904 slit lamp can provide high-quality examinations, even when you practice in more than one consultancy room, or perform field-based eye examinations.

The BA 904 can be used for hand-held operation to examine patients who cannot comfortably sit at a larger slit lamp, yet can be easily reassembled for more traditional joystick/headrest operation. Whichever method is used, the portability does not diminish the quality of the examination.

Both reliable and accurate, the BA 904 offers a Dioptic range of -8 to +8, an interpupillary range of 53 to 95. It offers a wider field of view, compared to other portable slit lamps available.

Both light and portable, the complete kit fits into the specially-designed light-weight suitcase. The integrated rechargeable NICAD battery allows up to 45 minutes of continuous use per battery (two batteries supplied).

## High-quality Haag-Streit optics

The BA 904 is the only portable slit lamp that incorporates World-class Haag-Streit optics. These superior optics offer high-resolution examinations; boasting 10x objectives, as standard and 16x objectives, as an option.

## Versatile

A dual-function slit lamp, the BA 904's ergonomic design enables both hand-held examinations and traditional-style examinations, using the easy-to-assemble head and chin rest.

## Portable

Extremely light and portable, the BA 904 effortlessly fits into a specially-designed suitcase, which is supplied as standard. It is easy to transport between consultancy rooms or field-based locations.

## Safe & comfortable

A well-balanced instrument, the BA 904 portable, hand-held slit lamp is both safe and comfortable for the patient. It filters out hazardous radiation, which means less damage to the retina when using a non-contact lens.

## Integrated battery

Its integrated, rechargeable NiCad battery allows up to 90 minutes examination time, away from the consultancy room. Two batteries are supplied, as standard.

## Two BA 904 models available

You can choose the portable slit lamp to suit your individual requirements from the two models that are available;

### Set A

BA 904 portable slit lamp, including; head and chin rest stand, two energy packs, charger, power supply and large case (Part Number: 5606011).

### Set B

BA 904 portable slit lamp, including; two energy packs, charger, power supply, parking unit and small case (Part Number: 5606013).

## Meets IEC 60601-1 3rd Edition & regulatory standards

The BA 904 has been updated to meet the stringent IEC 60601-1 3rd edition and regulatory standards. It has undergone rigorous safety, Electromagnetic Compatibility (EMC) and environmental testing, including; temperature range, vibration and humidity testing, as-well-as minor design improvements. All plastic parts are completely flame-retardant.



HIGH-QUALITY OPTICS



DUAL-FUNCTION SLIT LAMP



LIGHTWEIGHT & PORTABLE

## HAAG-STREIT UK

Edinburgh Way  
Harlow  
Essex  
CM20 2TT  
Phone (01279) 456261  
info@haag-streit-uk.com  
www.haagstreituk.com

## Members of HAAG-STREIT Group

HAAG-STREIT Holding AG  
www.haag-streit-holding.com

HAAG-STREIT Deutschland GmbH  
www.haag-streit.de

HAAG-STREIT AG  
www.haag-streit.com

IPRO GmbH  
www.ipro.com

SPECTROS AG  
www.spectros.ch

CLEMENT CLARKE Ltd.  
www.clement-clarke.com

HAAG-STREIT Medtech AG  
www.haag-streit-medtech.com

HAAG-STREIT UK  
www.haag-streit-uk.com

HAAG-STREIT France EURL  
www.haag-streit.fr

John Weiss Ltd.  
www.johnweiss.com

HAAG-STREIT Far East  
www.haag-streit-fareast.com

HAAG-STREIT USA  
www.haag-streit-usa.com

HAAG-STREIT Surgical GmbH  
www.haag-streit-surgical.com

Reliance Medical Inc.  
www.haag-streit-usa.com

Möller-Wedel GmbH & Co KG  
www.moeller-wedel.com

Asetronics AG  
www.asetronics.ch

Möller-Wedel Optical GmbH  
www.moeller-wedel-optical.com

ComLab AG  
www.comlab.ch