



WIZARD 210 & 220

Nasal Mask & Full Face Mask

Providing a comfortable fit is crucial for patient compliance, but it doesn't need to be complicated. The simple and effective design of WIZARD 210 & 220 allows the user to enjoy a relaxed night's sleep and encourages compliance to prescribed therapy. Its lightweight, soft cushion and ergonomic forehead pad allows for any pressure to be distributed evenly, improving the quality of therapy and providing maximum stability and comfort.



Better Mask, Better Experience

The simple but effective design of the WiZARD 210 & 220 brings users an easy good night's sleep. Its high quality silicon cushion and forehead pad provide a soft and smooth feel. To fit a wide range of facial contours, The WiZARD 210 is available in four comfortable sizes: XS, S, M and L; and the WiZARD 220 has six sizes: S, SW, M, MW, L, and LW (W type for wider nasal bridges). The WiZARD 210 & 220 features a lightweight and contoured design that ensures great comfort and an excellent seal, offering users a better experience for sleep therapy.



■ WiZARD 210
Nasal Mask

■ WiZARD 220
Full Face Mask



3 palpable fixation points for easy detachment/attachment of mask cushion.



Pressure or oxygen port with attached safety cap for pressure measurements or oxygen supplementation needs.



The model number and size is imprinted on frame and cushion for easy identification and replacement.



360° dual swivel elbow for freedom of movement during sleep.



Forehead pad with pressure distribution technology for excellent comfort.



Headgear clip for single handed operation for simple mask fitting and lifting.



Air-guiding vent channels for effective CO₂ wash-out and quiet air flow.



Headgear for comfortable skin climate.



Tubing quick release for convenient connection and disconnection of the air tubing.



Soft connection tubing with better flexibility for freedom of movement during sleep.

Ordering Information (For North America Only)

Order No.	Description	Order No.	Description	Accessories	
SM01008	WIZARD 210 Nasal Mask with Headgear, L	SM02011	WIZARD 220 Full Face Mask with Headgear, M	SM00001	WIZARD 210 Nasal Mask Cushion, L
SM01009	WIZARD 210 Nasal Mask with Headgear, M	SM02012	WIZARD 220 Full Face Mask with Headgear, S	SM00002	WIZARD 210 Nasal Mask Cushion, M
SM01010	WIZARD 210 Nasal Mask with Headgear, S	SM02016	WIZARD 220 Full Face Mask with Headgear, LW	SM00003	WIZARD 210 Nasal Mask Cushion, S
SM01011	WIZARD 210 Nasal Mask with Headgear, XS	SM02017	WIZARD 220 Full Face Mask with Headgear, MW	SM00030	WIZARD 210 Nasal Mask Cushion, XS
SM02010	WIZARD 220 Full Face Mask with Headgear, L	SM02018	WIZARD 220 Full Face Mask with Headgear, SW	SM00004	WIZARD 220 Full Face Mask Cushion, L

Accessories	Accessories	Accessories	
SM00005	WIZARD 220 Full Face Mask Cushion, M	SM00011	WIZARD 210/220 Forehead Pad, Standard
SM00006	WIZARD 220 Full Face Mask Cushion, S	SM00012	WIZARD 210/220 Forehead Pad, High
SM00025	WIZARD 220 Full Face Mask Cushion, LW	SM00028	WIZARD 210/220 Forehead Pad, Extra High
SM00026	WIZARD 220 Full Face Mask Cushion, MW	SM00016	Small Silicon Tubing
SM00027	WIZARD 220 Full Face Mask Cushion, SW	SM00034	WIZARD 210 Elbow Connector with Port Cap
		SM00035	WIZARD 220 Elbow Connector with Port Cap (M/L/LW)
		SM00036	WIZARD 220 Elbow Connector with Port Cap (S/SW/MW)
		SM00009	WIZARD 210/220 Headgear with Buckles
		SM00010	WIZARD 210/220 Buckles

HEADQUARTERS

Apex Medical Corp.
No. 9, Min Sheng St., Tu-Cheng,
New Taipei City, 23679, Taiwan
T +886 2 2268 5568
F +886 2 2268 6525
marketing@apexmedicalcorp.com

PAN-ASIA

TAIWAN

Apex Medical Corp.
No. 9, Min Sheng St., Tu-Cheng,
New Taipei City, 23679, Taiwan
T +886 2 2268 5600 /0800 054 133
F +886 2 2268 6526
customer@apexmedicalcorp.com

CHINA

Apex Medical (Kunshan) Corp.
No. 1368, ZiZhu Road, Kunshan,
Jiang Su Province, China
T +86 512 8617 1660 Ext 1150
F +86 512 8617 8663
customer.cn@apexmedicalcorp.com

INDIA

Apex Medicalcorp India Pvt. Ltd.
303, Third Floor Block 4A DLF
Corporate Park Phase-III,
Gurgaon, Haryana
T +91 11 6654 4138
F +91 11 6654 4052
enquiries.india@apexmedicalcorp.com

THAILAND

Apex Medical Thailand
65/100 Chamnanphenjati Business Center Bldg.
11/F Rama 9 Road, Huaykwang
Bangkok 10320
T +662 948 4848
F +662 948 4170
sales.thai@apexmedicalcorp.com

MIDDLE EAST • AFRICA

29, Misr Liltameer Building,
Stree#10, AREA#8, Sheraton Building, Heliop-
olis Cairo, Egypt.
T +202 2 2671303
F +202 2 2678649
raafat.wahby@apexmedicalcorp.com

EUROPE

SPAIN

Apex Medical S.L.
Elcano 9, 6ª planta 48008
Bilbao. Vizcaya. Spain
T +34 94 4706408
F +34 94 4706409
info@apexmedical.es

UK

Apex Medical Limited
Unit 21, Optima Park, Thames Road,
Crayford, Kent DA1 40X
T +44 (0) 132 2520 560
F +44 (0) 132 2554 182
enquiry@apexmedicalcorp.co.uk

FRANCE

Apex Medical France SAS
8 Blvd. de l'Épervière,
49000 Écouflant,
Pays De La Loire, France
T +33 2 72 79 27 13
F +33 2 85 73 07 67
info.france@apexmedicalcorp.com

GERMANY

SLK Vertriebsgesellschaft mbH
Lindenhorster Straße 38-40 44147
Dortmund
T +49 (0)231-925360-175
F +49 (0)231-925360-29
info@slk-gmbh.de
www.slk-gmbh.de

AMERICA

USA

Apex Medical USA Corp.
927 Mariner Street,
Brea, CA 92821, USA
T +1 714 671 3818
F +1 714 671 9040
sales.usa@apexmedicalcorp.com

Distributed by