

Application aids

Wall dispensers and dosage aids

Additional equipment for medical disinfection and cleaning ease the daily handling and dosing of the utilized products. Proper dosage aids can influence dripping and therefore reduce unnecessary wastage. Dispensers and spray heads enable a targeted application and exact released amount and hence an optimised consumption. With flexibly mountable holders, disinfection can conveniently be reached and ready for use in nearly anywhere. Dr. Schumacher makes it possible to comply with strict hygiene guidelines with these coordinated components.



AK 500/1000 PLUS

wall dispenser with aluminium casing

- autoclavable Euro dispenser with aluminium casing
- stainless steel pump can easily be removed without disassembling
- short or long stainless steel lever
- for 500 ml or 1 l bottles
- dosing 1.5 ml per stroke, can be reduced in three steps

Single unit	Unit	Size	REF
AK 500/1000 PLUS (long lever)			
wall dispenser	1	500 ml	00-901-005-01
wall dispenser	1	1 l	00-901-010-01
AK 500/1000 PLUS (short lever)			
wall dispenser	1	500 ml	00-903-005-01
wall dispenser	1	1 l	00-903-010-01



SPE 500/1000

armbar wall dispenser in plastic design

- for 500 ml or 1 l bottles
- closed design with window
- very stable, long plastic lever
- dosing 1.8 ml per stroke, can be reduced in three steps

Single unit	Unit	Size	REF
SP 500/1000			
wall dispenser	1	500 ml	00-918-005
wall dispenser	1	1 l	00-918-010



SPE TOUCHLESS 500/1000

touchless wall dispenser

- for 500 ml or 1 l bottles
- closed design with window
- inclusive drop catcher
- dosing 1.8 ml per stroke, can be reduced in three steps

Single unit	Unit	Size	REF
SP TOUCHLESS 500/1000			
wall dispenser	1	500 ml	00-911-005TL
wall dispenser	1	1 l	00-911-010TL



DESIPOINT

pre-assembled, mobile column stand for wall dispensers

- hygienic stainless steel construction
- quick and easy mounting
- for combination with all commercial dispensers (dispenser not included)

Single unit	Unit	Size	REF
DESIPOINT	1	160 cm	00-902-DS



DESIPOINT PRO

pre-assembled, mobile column stand for wall dispensers

- integrated TFT screen
- quick and easy mounting
- for combination with all commercial dispensers (dispenser not included)

Single unit	Unit	Size	REF
DESIPOINT PRO	1	150 cm	00-902-DSD



Dosage pumps

for hand and skin cleansing/disinfection products

- 1 dosage pump with 30 ml/stroke, fits 2 l bottle, 5 l or 10 l container adapter available on demand for smaller dosing (20 ml, 10 ml)
- 2 dosage pump with 3 ml/stroke, fits 500/1 l bottles
- 3 dosage pump with 2 ml/stroke, fits 500/1 l bottles

Single unit	Unit	Size	REF
Dosage pump with 30 ml/stroke	1	2l/5l/10l	00-906-000-01
Dosage pump with 3 ml/stroke	1	500ml/1 l	00-906-005-02
Dosage pump with 2 ml/stroke	1	500ml/1 l	00-906-005EX



Bottle holder

for hand and skin cleansing/disinfection products

- for wall mounting
- stable plastic construction
- suitable for 500ml/1 l dosage bottles

Single unit	Unit	Size	REF
wall holder	1	500ml/1 l	00-902-WHEX



DESCOFLEX bottle holder MH

for hand and skin cleansing/disinfection products

- to be put on bedsteads etc.
- made of chrome-plated steel
- easy fixing system
- requires no mounting materials
- suitable for 500 ml/1 l bottles

Single unit	Unit	Size	REF
bottle holder MH	1	500ml/1 l	00-902-MH



Dosage aids

for surface and instrument disinfectants and cleaners

- 1 dosage bottle 1 l with 60 ml measuring chamber
- 2 container opener
- 3 tap for 2 l bottle and 5, 10 and 30 l container

Single unit	Unit	Size	REF
dosage bottle	1	1 l with 60 ml measuring chamber	00-904-010
container opener	1	-	00-909-050
tap	1	2 l / 5 l / 10 l	00-907-000



Spray heads

for surface disinfectants and cleaners

- 1 spray head with adapter for 500 ml/1 l bottles
- 2 foam spray head with adapter for 500 ml/1 l bottles

Single unit	Unit	Size	REF
spray head	1	500 ml - 1 l	00-905-010
foam spray head	1	500 ml - 1 l	00-912-010



Dosage aids

for surface and instrument disinfectants and cleaners

- 1 measuring jug 3 l
- 2 measuring cup 40 ml

Single unit	Unit	Size	REF
measuring jug	1	3 l	00-906-030
measuring cup	1	40 ml	00-906-0004



DESCOMAT KOMPAKT DIGITAL

computer controlled dosage machine for disinfectants and cleaners

- especially user-friendly handling
- closed modular construction
- computerized with high functionality

DESCOMAT KOMPAKT DIGITAL is an innovative advancement of DESCOMAT KOMPAKT which established itself on the market during the last years. It is completely computer controlled and therefore easy to handle. Its outstanding technical functionality sets new standards for electronic dosage machines.

Especially user-friendly

- system status is shown on a four-row display
- easy to handle and monitor
- dosage rates: 0.25 – 5 % (optional up to 8 %), variable in 0.25 % steps
- pre-selection of use dilution amount continuously variable in 0.4 l steps
- two different application concentrations and issue quantities can be set up for one product (e.g. one for routine disinfection acc. to VAH, one for the case of an epidemic acc. to RKI)
- closed system construction that prevents bacterial contamination
- conforms to the guidelines of BAM (German Institute for material testing) and RKI (Robert Koch Institute)

Outstanding technical functionality

- easy maintenance
- USB-port for an easy upload of all usage-/failure data via PDA or notebook (optional)
- all functions of sensors and actuators displayed via LED (helpful for local technician)
- electronics covering made of acrylic glass
- lockable door, can be opened without tools
- made of stainless steel and other corrosion-resistant materials

Types and dimensions

- available as double-unit for two containers, allows for the dosing of two products at the same time
- available with green display (lightly illuminated) for rooms with daylight, alternatively with illuminated blue display e.g. for windowless rooms
- lightweight (8 kg), allows for an installation also on less sustainable walls
- dimensions (HWD): 30 x 20 x 17 cm

Single unit	Unit	Content	Art. Nr.
DESCOMAT Kompakt Digital	1	-	00-904D-KPE
DESCOMAT Kompakt double unit Digital	1	-	00-904D-KPD

Application aids overview

Application aids			
Wall dispenser	Unit	Size	REF
AK 500/1000 PLUS (long lever)	1	500 ml	00-901-005-01
	1	1 l	00-901-010-01
AK 500/1000 PLUS (short lever)	1	500 ml	00-903-005-01
	1	1 l	00-903-010-01
SPE 500/1000	1	500 ml	00-918-005
	1	1 l	00-918-010
SPE Touchless 500/1000	1	500 ml	00-911-005TL
	1	1 l	00-911-010TL
Wall dispenser accessories	Unit	Size	REF
Drop catcher AK 500 Plus	1	500 ml	00-901-205
Drop catcher AK 1000 Plus	1	1 l	00-901-210
DESIPOINT	Unit	Size	REF
DESIPOINT	1	160 cm	00-902-DS
DESIPOINT PRO	1	150 cm	00-902-DSD
Box and wall holders	Unit	Size	REF
Descoflex box holder DH	1	-	00-902-DH
Descoflex wall bottle holder WH	1	500 ml/1 l	00-902-WH
Descoflex mobile bottle holder MH	1	500 ml/1 l	00-902-MH
DESCO WIPES DT + ECO WIPES wall holder	1	-	00-902-EIM
DESCO WIPES DT SMALL wall holder	1	-	00-902-EIMK
one system PLUS - wall holder	1	-	00-902-OSP
one system basic - wall holder	1	-	00-902-OSB
Dosage aids	Unit	Size	REF
Dosage bottle	10	1 l	00-904-010
Container opener	10	-	00-909-050
Dosage pump with 30 ml/stroke	1	2 l/5 l/10 l	00-906-000-01
Tap	10	2 l/5 l/10 l	00-907-000
Dosage pump	1	500 ml/1 l	00-906-005-02
Spray head	10	500 ml/1 l	00-905-010
Foam spray head	10	500 ml/1 l	00-912-010
Measuring jug	1	3 l	00-906-030
Measuring cup	1	40 ml	00-906-0004
Disinfection tubs	Unit	Size	REF
Drill box	1	100 ml	00-908-001
1 liter	1	1 l	00-908-010
3 liter	1	3 l	00-908-030
4 liter	1	4 l	00-908-040
5 liter	1	5 l	00-908-050
10 liter	1	10 l	00-908-100
Dosage machines	Unit	Size	REF
DESCOMAT Kompakt Digital	1	-	00-904D-KPE
DESCOMAT Kompakt double unit Digital	1	-	00-904D-KPD
Dispenser buckets 5 liter	Unit	Size	REF
DESCO WIPES DT dispenser bucket empty	1	5 l	00-915-E001
DESCO WIPES DT dispenser bucket empty	6	5 l	00-915-E002
ECO WIPES dispenser bucket empty	1	5 l	00-915-EEW001
ECO WIPES dispenser bucket empty	6	5 l	00-915-EEW002
Dispenser buckets 3 liter	Unit	Size	REF
DESCO WIPES DT SMALL dispenser bucket empty	1	3 l	00-915-EK004
DESCO WIPES DT SMALL dispenser bucket empty	6	3 l	00-915-EK005