

**MORE PRODUCTS
MORE POSSIBILITIES
MORE PROFIT**

**CREATING
KIMAL
SPECIALITY PACKS**

KIMAL

BUILDING A KIMAL PROCEDURE PACK

Provides greater choice, visibility, efficiency and cost savings



Create it - supported by our specialist clinical knowledge



Plan it - with clear visibility of when you'll receive it




Timeline

- Start** - By observing your procedure we will create a bespoke pack for you
- Week 1** - Review and amend your suggested pack utilising the Pack+ app
- Week 2** - Non-sterile sample available for review
- Week 4** - Final sterile sample to approve
- Week 5** - Pack goes into production
- Complete** - Pack available for delivery after 6 weeks

Deliver it - on time with no unexpected surprises



Make it - with high quality products to meet your needs



Delivered with the Kimal Commitment

CONTENTS

WHY KIMAL?	4	SYRINGES	20
SPECIALITY PROCEDURE PACKS	5	NEEDLES	22
GETTING THE BASICS RIGHT	6	BOWLS	24
TROLLEY COVERS	8	SKIN PREPARATION	25
GOWNS AND TOWELS	10	SWABS	26
CUSTOMISED DRAPES	12	LABELS	27
MANIFOLDS AND LINES	14	TRAYS	28
WASTE MANAGEMENT	16	EQUIPMENT COVERS	29
SHARPS DISPOSAL	18	NOTES	30

WHY KIMAL?

Our mission is to be a;

**SPECIALITY PROCEDURE PACK MANUFACTURER
DRIVEN BY OUR CUSTOMERS' NEEDS.**

We are not a generic pack producer.

WHAT DOES THIS MEAN?

That we should differentiate our product offering by creating value added customised solutions tailored to our customers' needs.

SPECIALITY PROCEDURE PACKS

A Pack+ procedure pack is designed to provide the single use and consumable products for a specific procedure in a single package, the contents of which are specified by the user and can be hospital or clinic specific.

With over 3,000 components to choose from we have many specialised products to create a truly customised procedural solution.

When creating a pack for a customer it's important that we use our expertise to build a solution with our customers and they don't simply tell us what to include. It is essential, therefore, that you attend the procedure with the customer in order to fully understand their needs.

We have a core range of superior product solutions that provide opportunity to upspec existing or competitor generic packs.




IT'S IMPORTANT TO GET THE BASICS RIGHT

Whilst we have over 3,000 components to choose from it is important that you select from commonly used or cost-effective items to ensure that you create the optimum solution for your customer, whilst also ensuring we give yourselves the best possible opportunity to win new business.


Using the information in this booklet and information also available to you through Perform and the Pack+ App will ensure that you always select the correct commonly used items.

Once you have the basic and core components correct you can then focus on opportunities to Upspec and create Kimal speciality packs.



DELIVER ON TIME

Kimal recognise the importance of delivering on time, every time, in full. Through our 50 Years' experience we have developed core skills which have enabled us to work with our customers and actively manage the challenges of fluctuating demand through forecasted monthly usage.



Our **commitment**
to your **patient**

AWARD-WINNING CUSTOMER SERVICE

The customer is our number one priority and we are committed to consistently delivering high standards of customer care globally. It is this philosophy which has led to us to achieve the European Supply Chain Excellence Award for Customer Service in conjunction with an NHS Trust.



SPECIALIST CLINICAL KNOWLEDGE

Our team of Clinical Specialists and a customer focussed sales team are dedicated to understanding the needs of clinicians and will support you in improving efficiencies in your department.




KIMAL COMMITMENT

NO UNEXPECTED SURPRISES

We will work with you to deliver a customised procedure pack to your requirements. You will receive what you ask for and if anything changes, we will notify you first.

INNOVATIVE, HIGH QUALITY PRODUCTS TO MEET YOUR NEEDS

We work with our manufacturers to deliver quality products which will withstand the demands of the procedure.




TROLLEY COVERS

As an integral part of any pack the trolley cover provides a sterile field upon which to set up all other components, whilst also securing the contents of the pack itself.


DESCRIPTION	PRODUCT CODE	PREFERRED
140 X 210 BLUE	8360LONG/CN	YES
140 X 140	8360/CN	YES
140 X 175	8360BLUE/CN	YES
100 X 140	8368	YES
140 X 210 HEAVY DUTY	89009	
140 X 175 REINFORCED	83602	
140 X 280 REINFORCED	83603	
200 X 300 REINFORCED	V108-051-08	
140 X 210 REINFORCED	83601	YES
150 X 200 REINFORCED	S-51715	

TROLLEY COVERS

OPTIONS AVAILABLE



ZONE REINFORCED



FULLY REINFORCED

QUESTIONS TO ASK

- What size of trolley are you covering?
- What is the trolley used for other than holding the contents of the pack?
- Are there any heavy items likely to be placed on the trolley?
- Does the trolley cover need to absorb fluid?
- Do you require extra covers for additional trolleys?
- Do they require zoned or fully reinforced options?


GOWNS AND TOWELS

Gowns are an important part of your pack and vital to get right. Many Clinicians will have a preference of gown type and design, so it's important you understand the alternative and suitable options to recommend.

DESCRIPTION	PRODUCT CODE	PREFERRED
SMS LARGE	SGL	YES
SMS XL	SGXL	YES
SMS XL	8382BLUE/ VELCRO/XLONG	
SPUNLACE LARGE, 135CM	S-350033	
SPUNLACE XL, 145CM	S-350043	
SPUNLACE LARGE, 130CM	8380	YES
SPUNLACE XLARGE, 145CM	8382	YES
SMS LARGE REINFORCED	SGLR	YES
SMS XLARGE REINFORCED	SGXLR	YES
IMPERVIOUS REINFORCED	84577	
COTTON BLUE 58 X 44	8422	YES
BLUE 50 X 50	8415BLUE	YES
PAPER HAND TOWEL 57 X 33	8404	YES
WHITE 50 X 50	8415	YES

GOWNS AND TOWELS

OPTIONS AVAILABLE



STANDARD GOWN



REINFORCED GOWN

QUESTIONS TO ASK

- What properties do you want from a gown?
- Is there a large amount of fluid produced during the procedure?
- Do you have a preference for the type of fabric?
- Are you aware that spunlace is more expensive than SMS?
- How many hand towels would you like?
- How would you like the gowns and towels to be packed?

CUSTOMISED DRAPES

The SurgiGuard® range comprises a selection of materials and configurations to choose from. The opportunity to fully customise a draping platform is available to customers who have a specific preference or requirement.

DESCRIPTION	PRODUCT CODE	PREFERRED
RADIAL 240 X 330 4F 0W 3PLY	8403301	YES
RADIAL 345 X 240 4F 2W 3PLY	8423451	YES
RADIAL 225 X 305 4F 0W	50202406	YES
RADIAL 240 X 345 2W 4F SMS	8423452/B	YES
RADIAL 240 X 345 4F 2W 2PLY	50205211	YES
RADIAL 256 X 345 4F 2W SMS	V502-033-14	
PACING 220X 330 2F 2W SMS	V502-039-08	
PACING SMS SET (200 X 270 1F)	P504-004	YES
PACING 270 X 350 1F SMS	V502-009-14	
FEMORAL 220 X 360 2W 2F SMS	8223602/B	YES
FEMORAL 220 X 360 2F 1W SMS	V502-029-07	
FEMORAL 200 X 330 2F 1W 3PLY	8213301	YES
FEMORAL 220 X 330 2F 2W 3PLY	8223301	YES
FEMORAL 250 X 390 2F 2W 3PLY	50101106	YES
FEMORAL 250 X 350 2F 2W 2PLY	V502-031-13	
FEMORAL 200 X 330 2F 2W SMS	8223302/B	YES
FEMORAL 200 X 330 2F 1W 3PLY	8213301	YES
FEMORAL 200 X 330 2F 0W 3PLY	8203301	YES
EP 200 X 370 3F 2PLY	V502-003-16	

CUSTOMISED DRAPES

WHICH FABRIC TO CHOOSE?



2 ply



3 ply

LAMINATE

- Items can be placed on drape without causing strikethrough
- Fully impervious reducing risk of strikethrough
- Ideal for lengthy or 'bloody' procedures
- Absorbent & stays wet
- Molnlycke® customers

SMS

- Repellent for strikethrough protection
- Cost-effective
- More breathable
- Drapable
- Ideal for short or minimally invasive procedures
- KC customers



QUESTIONS TO ASK

- What type and length of procedure is the drape for?
- Do you know the size you require? What length and width?
- What fenestrations do you need?
- Size and location? Do you require adhesive or incise?
- Do you require windows? If so which side or where?
- Does the drape need an absorbent area?

MANIFOLDS & LINES

The hub of the fluid management system, managing the flow of contrast, saline and waste and allowing access to vascular pressures. A manifold can be supplied on its own or as part of a kit including fluid/contrast/manometer pressure line.

DESCRIPTION	PRODUCT CODE	PREFERRED
MANIFOLD.2 PORT RIGHT ON HALF	10001/HALF	
MANIFOLD 2 PORT R OFF HALF 250 PSI	10003/HALF	
MANIFOLD.2 PORT L OFF HALF 250PSI	10004/HALF	
MANIFOLD 3PORT RIGHT ON HALF 250PSI	10007/HALF	
MANIFOLD 3 PORT R OFF HALF 250PSI	10009/HALF	
MANIFOLD 3V RH/OF 200 RH STR OFF	10009/NAM	
MANIFOLD 2 VALVE RIGHT PORT OFF	80038200/DE	YES
MANIFOLD 2 PORT ON	80038202/DE	YES
MANIFOLD 3 PORT	80038301/DE	YES
MANIFOLD 3 PORT ON	80038303/DE	YES
MANIFOLD 2 PORT HP RIGHT PORT OFF	80037200/DE	
MANIFOLD 3 PORT HP RIGHT PORT OFF	80037301/DE	
MANIFOLD 3 PORT HP RIGHT OFF HALF	10109/HALF	

MANIFOLD KITS

DESCRIPTION	PRODUCT CODE	PREFERRED
STANDARD 2 PORT MANIFOLD KIT - VENTED SPIKE	8/NN/656608725	YES
STANDARD 2 PORT MANIFOLD KIT - SQUEEZE	8/NN/656608724	YES
RE-USEABLE SPIKE - FOR USE WITH SQUEEZE	H965904000391	YES

MANIFOLDS & LINES

LINES & CONTRAST DELIVERY

DESCRIPTION	PRODUCT CODE	PREFERRED
100CM MANOMETER LINE PLAIN	100MANFC/NS	YES
100CM MANOMETER LINE RED	100MANFR/NS	
100CM MANOMETER LINE BLUE	100MANFB/NS	
VENTED SPIKE LINE - CONTRAST - 180CM	9865	YES
VENTED CHAMBER LINE - CONTRAST - 180CM	9860	YES
NON-VENTED CHAMBER LINE - CONTRAST - 180CM	9860	YES
NON-RETURN VALVE	3191	YES

QUESTIONS TO ASK

- What type of manifold do you currently use? 2 port, 3 port, left hand, right hand, on/off?
- What length of manometer/pressure/contrast line would you need?
- What length of manometer line is needed with which connection and is it required to be plain or blue/red stripe?
- What type of drip chamber and spike is needed, vented/non-vented or contrast/flush line?
- Which lines do you require to be attached to which port of the manifold?
- What length of high pressure catheter extension is needed?

WASTE MANAGEMENT

Depending on current practice a customer may select from a wide range of options to safely dispose of contaminated waste fluid during the procedure.

DESCRIPTION	PRODUCT CODE	PREFERRED
BAG DRAINAGE 2000ML	0972/A	YES
FLUID SAFE BAG (K62/FCB)	0972/A + 3191	YES
CONTRAST GIVING SET / NEEDLE FREE (K90/36)	0972/A + 0040	YES
KIMAL FLUID COLLECTION DEVICE, LIDDED BOWEL WITH SOLIDIFIER	8543/FCD	YES
NAMIC PROTECTION STATION	8/NN/64000934	YES
NAMIC PROTECTION STATION PLUS	8/NN/656206817	YES
NAMIC PROTECTION STATION WITH DRAPE CLIP	8/NN/656608747	YES

WASTE MANAGEMENT

OPTIONS AVAILABLE



**FLUID
COLLECTION
DEVICE**



CLOSED SYSTEMS

QUESTIONS TO ASK

- How do you dispose of fluid waste currently?
- How often do you see open bowls splash back or spill?
- What potential risks are you aware of that are posed by contaminated waste?
- What benefits do you think using a closed system would bring?

SHARPS DISPOSAL

All procedures involve the use of needles, scalpels or blades that need to be safely disposed of. Needle stick injury is a risk faced by all healthcare workers who come into contact with these items. The method of disposal in each hospital may differ, dependant upon local policy and protocol or just outdated practice.

DESCRIPTION	PRODUCT CODE	PREFERRED
NEEDLE COUNTER FOAM MAGNET RED	8623	YES
NEEDLE CONTAINMENT DEVICE, KIMAL SOLUTIONS	8538/NCD	YES
NEEDLE COUNTER 20-COUNT DOUBLE MAGNET RED	85069	
NEEDLE COUNTER 40-COUNTFOAM MAGNET YELLOW	84982	
NEEDLE COUNTER SHARPS CONTAINER BOX LARGE YELLOW	8543	
NEEDLE COUNTER SHARPS CONTAINER BOX MEDIUM YELLOW	8/ASPEN/3	YES
NEEDLE COUNTER SHARPS PAD PURPLE TURTLE	8519	YES
NEED COUNTER SHARPS PAD ADHESIVE AND FOAM, CARD	8525	YES

SHARPS DISPOSAL

OPTIONS AVAILABLE



**NEEDLE
CONTAINMENT
DEVICE**



**CLOSED
CONTAINERS**

QUESTIONS TO ASK

- How do you dispose of sharps currently?
- How many needles/blades do you use?
- How do you keep track of the number of needles used during a procedure?
- Have you seen needles get lost or dropped? How often does this occur?
- Has anyone you know experienced a needle stick injury? What are the risks you think this may carry?
- How is current practice safeguarding yourself or colleagues?

SYRINGES

Syringes are used for a variety of things during any procedure including drawing up fluid or drugs and administering contrast.

DESCRIPTION	PRODUCT CODE	PREFERRED
BD SYRINGE 1ML, LUER LOCK	6801/BD	
BD SYRINGE 1ML, LUER SLIP	6901/BD	
SYRINGE 2ML, LUER LOCK	6802/TER	
SYRINGE 2ML, LUER SLIP	6902	
SYRINGE 3ML, LUER LOCK	6803	YES
SYRINGE 5ML, LUER SLIP	6905	YES
BD SYRINGE 5ML, LUER SLIP	6905/BD	
SYRINGE 5ML, LUER LOCK	6805	YES
BD SYRINGE 5ML, LUER LOCK	6805/BD	
BD SYRINGE 10ML, LUER LOCK	6810/BD	YES
SYRINGE 10ML, LUER LOCK	6810	
BD SYRINGE 10ML, LUER SLIP	6910/BD	YES
SYRINGE 10ML, LUER SLIP	6910	YES
BD SYRINGE 20ML, LUER LOCK	6820/BD	YES
SYRINGE 20ML, LUER LOCK	6820	YES
BD SYRINGE 20ML, LUER SLIP	6920/BD	
SYRINGE 20ML, LUER SLIP	6920	YES
SYRINGE 30ML, LUER LOCK	6830	
SYRINGE 10ML PP, LUER LOCK, DRUG SAFETY BLUE	6771BLUE	
SYRINGE 10ML PP, LUER LOCK, DRUG SAFETY GREEN	6771GREEN	
SYRINGE 10ML PP, LUER LOCK, DRUG SAFETY RED	6771RED	

SYRINGES

DESCRIPTION	PRODUCT CODE	PREFERRED
SYRINGE 10ML PC, LUER LOCK PLUNGER BUE	80505103/A	YES
SYRINGE 12ML PC, LUER LOCK CONTROL 3 RINGS	80087007/A	YES
SYRINGE 10ML PC, LUER LOCK CONTROL 3 RINGS ROTATING ADAPTER	6575/NAM	YES
SYRINGE 10ML PC, CONTROL	6590	YES
SYRINGE 10ML PC, LUER LOCK WHITE	6610/DMC	YES
SYRINGE 10ML PC, LUER LOCK GREEN	6610/DMC/G	YES
SYRINGE 10ML PC, LUER LOCK RED	6610/DMC/R	YES
SYRINGE 10ML PC, LUER LOCK BLUE	6610/DMC/B	YES
SYRINGE 10ML PC, LUER LOCK CONTROL PLUNGER	80095107/DE	
SYRINGE 10ML PC, LUER LOCK CARDIO CONTROL	6590	

QUESTIONS TO ASK

- What aspect of the procedure do you need syringes for?
- Drawing up or administering drugs? If so what volume/connection?
- Flushing/preparing products? If so what volume/connection?
- Administering contrast?
- Which handle/grip would you prepare?
- Will the syringe connect to a manifold?
- Would coloured syringes be useful to denote differing contents?

NEEDLES

Needles are required either to draw up or administer medication or facilitate the placement of wires during a procedure. For each requirement the gauge, length and type may vary.

DESCRIPTION	PRODUCT CODE	PREFERRED
NEEDLE 14G X 8CM GREEN	6526	YES
NEEDLE 18G X 50MM HYPO PINK	6517	YES
NEEDLE 18G X 3.8CM HYPO PINK	6527	YES
NEEDLE 19G X 3.8CM HYPO BEIGE	6528	YES
NEEDLE 20G X 3.8CM HYPO YELLOW	6530	YES
NEEDLE 20G X 7CM HYPO	6518	YES
NEEDLE 21G X 3.8CM HYPO GREEN	6531	YES
NEEDLE 21G 5CM HYPO GREEN	6539	YES
NEEDLE 21G X 8CM GREEN	6525	YES
NEEDLE 22G X 32MM HYPO BLACK	6523	YES
NEEDLE 22G 4CM HYPO BLACK	6532	YES
NEEDLE 23G X 3CM HYPO BLUE	6533	YES
NEEDLE 25G X 5/8" HYPO ORANGE	6535	YES
NEEDLE 25G X 2.5CM HYPO ORANGE	6536	YES

INTRODUCER NEEDLES

DESCRIPTION	PRODUCT CODE	PREFERRED
NEEDLE 18G X 7CM SELDINGER	6178	
NEEDLE 18G X 7CM (70MM) INTRODUCER	6180	YES
NEEDLE 19G X 7CM SINGLE WALL	6185	
NEEDLE 18G X7CM SING WALL CORDIS A4 PINK	6280	YES
NEEDLE 19G X7CM SING WALL CORDIS A4	6281	
SELDINGER NEEDLE 18G X 70MM	6298/A	

NEEDLES

SAFETY & DRAWING UP NEEDLES

DESCRIPTION	PRODUCT CODE	PREFERRED
BD ECLIPSE SAFETY HYPODERMIC NEEDLE 18G X 1.5"	6689	
BD ECLIPSE SAFETY HYPODERMIC NEEDLE 25GX5/8"	6693	
BD ECLIPSE SAFETY HYPODERMIC NEEDLE 23G X 1.25"	66018	
MAGELLAN™ SAFETY HYPODERMIC NEEDLE 18G X 1 1/2"	66001	
MAGELLAN™ SAFETY HYPODERMIC NEEDLE 25G X 5/8"	66012	
MAGELLAN™ SAFETY HYPODERMIC NEEDLE 23G X 1"	66011	
DRAWING UP NEEDLE BLUNT	8851	YES
NEEDLE BD DRAWING UP	8850	
NEEDLE 19G X 3.8CM FILTER 5 MICRON	6557	
NEEDLE 18G X 4CM FILTER BLUNT BEVEL	6603	
NEEDLE 18G BLUNT FILL WITH 5 MICRON FILTER	66027	YES

QUESTIONS TO ASK

- What aspect of the procedure do you need needles for?
- Which gauge and length would you prefer?
- Do you currently or have you considered switching to safety needles?
- Are you using the introducer needle to puncture the radial or femoral artery?

BOWLS

Used to store fluid or components, the requirement for bowls may vary between procedures. The Kimal range contains all common sizes and colours required.

DESCRIPTION	PRODUCT CODE	PREFERRED
GALLIPOT 60ML	36519084	YES
120ML GALLIPOT	36519086	YES
500ML BLUE	8581	YES
500ML CLEAR	8581/CLEAR	YES
250ML BLUE	8498	YES
250ML CLEAR	36519085	YES
1000ML BLUE	8499	YES
1500ML BLUE	8583	YES
6000ML BLUE	8504	YES
GUIDEWIRE BOWL, KIMAL SOLUTIONS	8606	YES
GUIDEWIRE BOWL 30CM RECTANGLE	8607/DE	

QUESTIONS TO ASK

- Which fluids do you usually decant into separate bowls and how much?
- How do you distinguish between the contents of different bowls?
- How do you currently prepare and store wires or catheters during the procedure?

SKIN PREPARATION

The skin is cleaned prior to any procedure. The method of cleaning will depend on the procedure or customer protocol and will vary from the use of simple sponges to prefilled applicators.

DESCRIPTION	PRODUCT CODE	PREFERRED
SPONGE STICK 20CM YELLOW	8555	YES
FOAM CUBE PREP SPONGE 5 X 4	84581	
SPONGE STICK LARGE 20CM	8611	
SPONGE STICK MEDIUM 15CM	8554	YES

QUESTIONS TO ASK

- How do you currently prepare the skin?
- Which solution do you use?
In which container?
- Do you prefer sponge sticks or individual sponges?
- Do you require a sponge holder?
- How many would you usually use?

SWABS

Swabs are commonly used to clean or dry patient skin. The size, quantity and colour may vary depending on the use in each procedure.

DESCRIPTION	PRODUCT CODE	PREFERRED
SWAB GAUZE 10 X 10 12PLY BLUE	8353BLUE/1	
SWAB GAUZE 10 X 10 12PLY TIED5	8320/5	YES
SWAB GAUZE 10 X 10 32PLY	8314/5	
SWAB GAUZE 5 X 5 4PLY 5PACKS	8336/5	YES
SWAB GAUZE 7.5 X 7.5 12PLY	8353/DE	
SWAB GAUZE 10 X 10 XR 12PLY 5PACK	8325/5	YES
SWAB GAUZE 30 X 30 XR 12PLY 5PACK	8332/5	

QUESTIONS TO ASK

- What do you use the swabs for?
- What size and how many would you need?
- How do you keep track of the number of swabs used?
- Does the procedure involve any surgical incision?

LABELS

Commonly used to label syringes or gallipots, labels help easily identify the medication or solution each contains. Labels are available as a sheet or pre-applied to syringes.

DESCRIPTION	PRODUCT CODE	PREFERRED
LABEL WATERPROOF 45 X 10 PURPLE	1121PURPLE	
LABEL WATERPROOF 45 X 10 BLUE	1121BLUE	
LABEL WATERPROOF 45 X 10 RED	1121RED	
LABEL WATERPROOF 45 X 10 YELLOW	1121YELLOW	
LABEL WATERPROOF 45 X 10 GREEN	1121GREEN	
LABEL SET WATERPROOF	1160/NEW	YES
LABEL WATERPROOF 45 X 10 BROWN	1121BROWN	
LABEL WATERPROOF 45 X 10 WHITE	1121WATER	

QUESTIONS TO ASK

- How do you currently identify the contents of syringes if using more than one?
- What medication/fluid will you use during the procedure?

TRAYS

A tray is generally used to form the base of a procedure pack and provide stability. Some trays may also be selected to specifically hold or protect products during the procedure or whilst in transit.

DESCRIPTION	PRODUCT CODE	PREFERRED
TRAY FULL DEEP 28 X 25	6019	YES
TRAY HALF DEEP 24.5 X 13.5	6018	YES
TRAY 17X 5 X 13.5 CLEAR	6015	YES
TRAY HALF SHALLOW 24.5 X 13.5	6017	YES
TRAY DEEP 42 X 30	6023	

QUESTIONS TO ASK

- Do you currently use a tray to hold any products during the procedure?

EQUIPMENT COVERS

During the procedure to prevent cross contamination or de-sterilisation a number of pieces of equipment will be covered. The size and requirement for covers will vary depending on the equipment used in each lab.

DESCRIPTION	PRODUCT CODE	PREFERRED
COVER FLUORO SMALL RED	8099	YES
COVER FLUORO MEDIUM BLUE	8100	YES
COVER FLUORO LARGE GREEN	8101	YES
COVER 102 X 102CM	8111	YES
COVER EQUIPMENT 76 (OPENING) X 91CM	8252	YES
COVER 100 X 75CM	8257	YES

QUESTIONS TO ASK

- What equipment do you need to cover during the procedure?
- Do you think your equipment adequately protected?

NOTES



KIMAL