



varimed[®]
Medical Furniture
for Examination and Treatment Areas



varimed®

Couches, stools, medi-net tables, examination seats –
of indispensable help for the daily routine

Long-wearing upholstery with
leatherette cover.

Rounded, reinforced
edges, seamless
upholstery.

All examination couches with spacer
block between frame and upholstery -
no risk of crushing fingers!

Levelling feet to
compensate unevenness
of floors.

Easy and fast height
adjustment by
"anatomic-Lift"
gas spring.

Friction-braked double,
plastic castors resistant
to abrasions.

Colours:
upholstery mandarin 98
powder coating grey aluminium RAL 9007



2 ratchets completely made of metal - for longevity.

Tubular steel frame of solid construction, 40 x 20 mm, feet of couches, 30 x 30 mm.

Drawers of medi-net trolley move on ball-bearing precision rails for an easy running; designed also for high working stress.

Tubular steel frame of solid construction, 20 x 20 mm.

High-quality powder coating of extreme knock and scratch resistance.

Double, plastic castors resistant to abrasions.



varimed® Couches

Examination couches
211.1180.0 and 211.1190.0



Colours:
211.1180.0
upholstery topaz 105, powder coating grey white RAL 9002

211.1190.0
upholstery lime 99, frame chromium-plated

Examination couch 211.1080.0



Colours:
upholstery mandarin 98
powder coating grey aluminium RAL 9007

ECG examination couch 211.1310.0









Colours:
upholstery opal green 100
powder coating white aluminium RAL 9006

Bariatric couch 211.1330.0



Colours:
upholstery caribbean blue 82
frame chromium-plated

All examination couches with spacer block between frame and upholstery!

Model No.	Description
211.1080.0 	Examination couch €€ rectangular steel tube frame, powder coated RAL 9002, RAL 9006 or RAL 9007, 2-section 50 mm thick foam rubber upholstery, levelling feet, head section adjustable up to +40°. Safe working load: 225 kg.
211.1180.0 	Examination couch €€ rectangular steel tube frame, powder coated RAL 9002, RAL 9006 or RAL 9007, 2-section 70 mm thick foam rubber upholstery, levelling feet, head section adjustable up to +40°. Safe working load: 225 kg.
211.1190.0	Examination couch €€ chromium-plated frame, otherwise like 211.1180.0
200.0970.0 200.1030.0 200.6500.0 	Accessories: Paper roll holder for paper roll up to 500 mm width Paper roll , 500 mm wide Surcharge for coloured frame acc. to colour chart, for 211.1080.0 and 211.1180.0 Surcharge for head section with nose inlet (for 211.1180.0 and 211.1190.0) Filling for nose inlet
211.1310.0 	ECG examination couch €€ 800 mm wide, rectangular steel tube frame, powder coated RAL 9002, RAL 9006 or RAL 9007, with lateral bracing, 70 mm thick foam rubber upholstery with 2-section cover, levelling feet, head section adjustable up to +40°. Safe working load: 225 kg.
211.1320.0	ECG examination couch €€ chromium-plated frame, otherwise like 211.1310.0
200.0940.0 200.1030.0 200.5100.0 200.5110.0 200.5120.0 200.6500.0 	Accessories: Paper roll holder for 800 mm wide couches Paper roll , 500 mm wide Longitudinal bracing , powder coated RAL 9002, RAL 9006 or RAL 9007, for 211.1310.0 Longitudinal bracing , chromium-coated, for 211.1320.0 Longitudinal bracing , powder coated according to colour chart Surcharge for coloured frame acc. to colour chart, for 211.1310.0
211.1330.0 	Bariatric couch €€ 800 mm wide, chromium-plated rectangular steel tube frame, with transverse and longitudinal bracing, 2-section mattress, 70 mm thick foam rubber upholstery, levelling feet, head section adjustable up to +40° Safe working load: 300 kg. Mobile version and special widths/lengths not available!
200.0940.0 200.1030.0	Accessories: Paper roll holder for 800 mm wide couches Paper roll , 500 mm wide
200.0000.0 200.0010.0 200.0020.0 200.0030.0	Surcharge for couches 211.1180.0 – 211.1320.0 with special design Different height up to max. 830 mm (from 700 mm height upwards, including transverse and longitudinal bracing) Mobility: 4 swivel-mounted castors diam. 125 mm, 2 castors with total locking device, antistatic, 1 castor with directional locking, with revolving bumpers diam. 150 mm, longitudinal bracing. All castors are resistant to abrasions. (please indicate desired final height of couch!) Different width of +/- 150 mm (max. width of 900 mm) Different length of +/- 150 mm

If both special height and width are required, surcharge will only be applicable once.

Note:
 For orthopaedic and massage treatment, couches with transverse and longitudinal bracings are recommended!

4 colours for leatherette upholstery in "classic" quality without surcharge.
18 colours for leatherette upholstery in "comfort" quality (200.8510.0/200.8520.0).
 Please see pages 18 and 19.

Standard powder-coating colours for frame: RAL 9002, RAL 9006 or RAL 9007.
 Please indicate desired colour upon order placement.

varimed® Couches, mobile Recovery couches, mobile

Examination & massage couch 211.9351.0



Colours:
upholstery blossom white 84
frame chromium-plated

Recovery couch 211.9501.0
1 Storage basket 200.5500.0




Colours:
upholstery caribbean blue 82
frame chromium-plated

All examination couches with spacer block between frame and upholstery!

Model No.

Description


211.9301.0



Examination & massage couch €€

with transverse and longitudinal bracing, rectangular steel tube frame, powder coated RAL 9002, RAL 9006 or RAL 9007, with 4 swivel-mounted castors diam. 125 mm, 2 castors with total locking device, antistatic, 1 castor with directional locking, with revolving bumpers diam. 150 mm, 70 mm foam rubber upholstery with 2-section cover, head section adjustable up to +40°. Safe working load: 225 kg.


211.9351.0



Examination & massage couch €€

chromium-plated frame, otherwise like 211.9301.0.


211.9501.0



Recovery couch €€

chromium-plated rectangular steel tube frame, with transverse and longitudinal bracing, with 4 swivel castors diam. 125 mm, 2 castors with total locking device, antistatic, 1 castor with directional locking, with revolving bumpers diam. 150 mm. 70 mm thick foam rubber upholstery with 2-section cover, head section adjustable up to +40°. Equipped with 2 side guards with attachments, chromium-plated, can be lowered down to height of upholstery. Safe working load: 225 kg.

211.9511.0



Recovery couch €€

chromium-plated frame, 680 mm height, otherwise like 211.9501.0.

Accessories:

for couches 211.9301.0 – 211.9511.0

200.0970.0



Paper roll holder

for paper roll up to 500 mm width.


200.1030.0



Paper roll

500 mm wide.


200.4900.0



1 Side guard and 1 side rail

side guard laterally adjustable on side rail, chromium-plated, with clamp, fitted to rectangular steel tube frame couches, may be lowered down to height of upholstery.

200.4950.0




1 Side guard with clamp

screwed, rigid connection, fitted to rectangular steel tube frame couches, chromium-plated, may be lowered down to height of upholstery.

200.0020.0

Surcharge for 800 mm width of couch (211.9301.0 to 211.9511.0)

200.5500.0



Storage basket

formed to nest onto transverse or longitudinal bracing, basket 520 x 490 x 140 mm, chromium-plated.

200.6500.0

Surcharge for coloured frame

acc. to colour chart, for 211.9301.0.

All castors are non-marking and resistant to abrasions.

4 colours for leatherette upholstery in "classic" quality without surcharge. 18 colours for leatherette upholstery in "comfort" quality (200.8520.0). Please see pages 18 and 19.

Standard powder-coating colours for frame: RAL 9002, RAL 9006 or RAL 9007. Please indicate desired colour upon order placement.



Please refer to page 4 of our catalogue no. 81 for further information on this product.

varimed[®] Height adjustable Examination and Treatment couch



- **Large and ergonomic height adjustment range**
- **No displacement of table-top during height adjustment**
- **Moulded base provides generous legspace**




Description

Examination and treatment couch € 211.4100.0

hydraulic height adjustment 500 to 920 mm by means of foot lever.
 Servo-assisted head section adjustment from -20° to $+40^{\circ}$.
 2-section upholstery 70mm thick, leatherette cover acc. to SCHMITZ colour chart.
 Powder coated chassis in white aluminium RAL 9006, grey aluminium RAL 9007 or grey white RAL 9002. Silver-grey plastic chassis cover.
 With 4 levelling feet. Safe working load: 225 kg

Examination and treatment couch € 211.4110.0

800 mm wide, otherwise like 211.4100.0.



Examination and treatment couch € 211.4200.0

electromotive height adjustment from 500 to 920 mm by means of foot control unit, with safety ON/OFF switch to turn control element on or off, 2.40 m connecting cable with mains plug, otherwise like 211.4100.0.
 Upon request, hand control instead of foot control unit.

Examination and treatment couch € 211.4210.0

800 mm wide, otherwise like 211.4200.0.



1

Safety ON/OFF switch

turns control element on or off.

Accessories:

2

Mobility 200.0250.0 4 closed double castors \varnothing 100 mm, 2 of which el. conductive, central lock per axis, 1 castor with directional lock. (Height adjustment range remains the same.)

Mobility, retrofit kit 211.9940.0

for examination and treatment couches with a width of 650 mm. (w/o ill.)



Mobility, retrofit kit 211.9940.1

for examination and treatment couches with a width of 800 mm. (w/o ill.)



3

Side guards, pair 200.0400.0

length: 990 mm, folding down to below top of mattress, powder coated white aluminium RAL 9006.

4

Holder for foot control 200.0350.0 for Examination and Treatment couches 211.4200.0 and 211.4210.0, of steel, powder coated white aluminium RAL 9006, for mounting at the long side of the couch (retrofitable)

Paper roll holder 200.0930.0 for paper roll up to 590 mm wide.

Paper roll 200.1030.0 500 mm wide.

Edge protection 200.0410.0 of plastic, grey, at all four upholstery corners (w/o ill.)

Head section with nose outlet 200.0500.0 (w/o ill.)

Filler for nose outlet 200.0510.0 (w/o ill.)

Side rail (piece) 200.4402.0 for Examination and Treatment couches 211.4100.0, 211.4110.0, 211.4200.0 and 211.4210.0, mounting either side, 25 x 10 x 750 mm, max. charge 30 kg, cannot be combined with side guards 200.0400.0. (w/o ill.)

ECG shield 200.0251.0 with one potential equalization sockets each LHS and RHS as well as screening plate to reduce disturbance, suitable for 650 mm wide couches. (w/o ill.)



ECG shield 200.0252.0 for 800 mm wide couches, otherwise like 200.0251.0. (w/o ill.)



Surcharge for coloured chassis 200.6530.0 acc. to colour chart.

4 colours for leatherette upholstery in "classic" quality without surcharge. 18 colours for leatherette upholstery in "comfort" quality (200.8520.0). Please see pages 18 and 19. Standard powder-coating colours for frame: RAL 9002, RAL 9006 or RAL 9007. Please indicate desired colour upon order placement.



varimed® Universal Examination and Treatment couches



Description

Universal examination couch € 211.4300.0

electromotive height adjustment from 500 to 920 mm by means of foot control unit.
With safety ON/OFF switch to turn control element on or off. Mid-section can be used as back rest, manual adjustment of head, back and seat section by 35° each, stainless steel rinsing basin, clamps for leg supports. 70 mm thick foam rubber upholstery, leatherette cover acc. to SCHMITZ colour chart. Powder coated chassis in white aluminium RAL 9006, grey aluminium RAL 9007 or grey white RAL 9002. Silver-grey plastic chassis cover. With 4 levelling feet.
2.40 m connecting cable with mains plug. Upon request, hand control instead of foot control unit.
Safe working load:
in horizontal position for general examinations: 225 kg
in gyn. position: 135 kg
Mobile version available as an option (model no. 200.0250.0).
Special dimensions not possible.

1 **Safety ON/OFF switch** turns control element on or off.

2

Paper roll holder 200.0930.0

for 211.4300.0, for paper rolls up to 590 mm width, fitting at head section.

Paper roll 200.1030.0

500 mm wide, 50 m roll. (w/o ill.)

3

Side rail (piece) 200.4401.0

at head section, 25 x 10 x 480 mm.



Description



4

Infusion pole 200.4600.0

stainless steel, upper part with 4 hooks, attachment clamp and ball-and-socket joint adjustment, range of adjustment 410 mm, fitting to side rail 200.4401.0 (not included), safe working load: 2 kg (per hook).



5

Mobility 200.0250.0

4 closed double castors Ø 100 mm, 2 of which el. conductive, central lock per axis, 1 castor with directional lock. (Height adjustment range remains the same.)

Mobility, retrofit kit 211.9940.0

for examination and treatment couches with a width of 650 mm. (w/o ill.)



6

Leg support (piece) 202.2270.0

Goepel type, silver-grey integral foam pad, with fastening strap, stainless steel 18/10 rod.

7

Arm rest 202.1070.0

plastic board, integral foam pad, antistatic, with attachment strap, stainless steel clamp, one-touch swivelling through 180° and horizontal as well as vertical adjustment, height adjustable, width 150 mm, length 600 mm. Mounting to side rail 200.4401.0 (not included).



7

8

Holder for foot control 200.0350.0

of steel, powder coated white aluminium RAL 9006, for mounting at the long side of the couch (retrofittable)

Surcharge for coloured chassis 200.6540.0

acc. to colour chart, for 211.4300.0.



8

4 colours for leatherette upholstery in „classic“ quality without surcharge. 18 colours for leatherette upholstery in „comfort“ quality (200.8520.0). Please see pages 18 and 19.

Standard powder-coating colours for frame: RAL 9002, RAL 9006 or RAL 9007. Please indicate desired colour upon order placement.





varimed®

Swivel stools, Examination seats, Medical furniture

Colours:
222.2760.0 upholstery milky coffee 102
222.2810.0 upholstery raspberry red 103



Colours:
upholstery caribbean blue 82

Model No.	Description
222.2810.0 	Swivel stool 5-arm plastic base, black, chromium-plated column, friction braked double castors for smooth/hard floors, height adjustment by means of "anatomic-lift" gas spring, with hand ring release, chromium-plated column, leatherette upholstered seat. Please specify double castors for use on carpeting upon order, if required. Safe working load: 135 kg.
200.8000.0	Surcharge for powder coated base according to colour chart; chromium-plated column.
222.2760.0 	Swivel stool 5-arm polished aluminium base, chromium-plated column, friction-braked double castors for smooth/hard floors, height adjustment by means of "anatomic-lift" gas spring, with hand ring release, leatherette upholstered seat. Please specify double castors for use on carpeting upon order, if required. Safe working load: 135 kg.
222.1550.0 	Examination seat specially designed upholstered seat with leatherette cover acc. to colour chart, height adjustment from 490 to 690 mm by means of gas spring, with foot ring release. 5-arm base powder coated in grey aluminium RAL 9007, with antistatic plastic castors Ø 50 mm for smooth/hard floors, with perimetral black bumper protection. Safe working load: 135 kg.
222.1570.0 	Examination seat specially designed upholstered seat with leatherette cover acc. to colour chart, height adjustment from 490 to 690 mm by means of gas spring, with foot ring release, foam rubber backrest with leatherette cover, 5-arm base powder coated in grey aluminium RAL 9007, with antistatic plastic castors Ø 50 mm for smooth/hard floors, with perimetral black bumper protection. Safe working load: 135 kg.
200.0011.0	Surcharge for brake loaded castors braked when loaded, released when you come off the stool for Examination seats 222.1550.0 and 222.1570.0

4 colours for leatherette upholstery in "classic" quality without surcharge.
18 colours for leatherette upholstery in "comfort" quality (surcharge 200.8530.0 for swivel stools 222.2810.0 and 222.2760.0/200.8540.0 for examination seat 222.1550.0/200.8550.0 for examination seat 222.1570.0).
Please see pages 18 and 19.

221.1860.0



Upholstery "classic" silver-grey 55

Model No.

Description

221.1860.0

Examination/Blood sampling chair €€


base of rectangular tube, powder coated in RAL 9002, RAL 9006 or RAL 9007, seat, backrest and armrests upholstered in "classic" black 57 leatherette. Safe working load: 225 kg.

Surcharge special coloured finish of frame acc. to colour chart.

**3 more colours for leatherette upholstery in "classic" quality (surcharge 200.8500.0).
18 colours for leatherette upholstery in "comfort" quality (surcharge 200.8560.0).
Please see pages 18 and 19.**

**Standard powder-coating colours for frame: RAL 9002, RAL 9006 or RAL 9007.
Please indicate desired colour upon order placement.**


200.6610.0



242.2350.0



242.2350.0



Drip-feed stand €€/93/42/EEC

18/10 stainless steel frame, disassembled, 5-arm base, upper section with 4 hooks, drip glass, double plastic castors diam. 50 mm, antistatic, rust-proof; height adjustment from 1,300 to 2,120 mm, one-hand action. Safe working load: 2 kg (per hook).

varimed[®] medi-net General purpose tables Equipment tables

medi-net equipment table 233.6060.0,
 frame grey aluminium RAL 9007



medi-net general purpose table 233.6010.0
 medi-net general purpose table 233.7010.0,
 frame grey white RAL 9002



Description


medi-net equipment table**

rectangular steel tube frame, collapsible, powder coated grey white RAL 9002, white aluminium RAL 9006 or grey aluminium RAL 9007, drawer with pre-assembled black galvanized ball bearing precision runners, table top, shelf and drawer front in pearl white 60 plastic veneers, handle for drawer front powder coated in grey white RAL 9002, white aluminium RAL 9006 or grey aluminium RAL 9007, 4 double castors diam. 75 mm.

The equipment table 233.6060.0 cannot be fitted with additional drawers.

Model No.

233.6060.0



Frame powder coated,
 1 drawer, height 750 mm

Description

medi-net general purpose tables**

rectangular steel tube frame, collapsible, chromium-plated or powder coated grey white RAL 9002, white aluminium RAL 9006 or grey aluminium RAL 9007, drawer with pre-assembled black galvanized ball bearing precision runners, table top, lower shelf and drawer front in pearl white 60 plastic veneers, handle for drawer front powder coated in grey white RAL 9002, white aluminium RAL 9006 or grey aluminium RAL 9007, 4 double castors diam. 75 mm.

Fitting of 2 drawers max. in addition to the standard drawer for table models 233.6010.0/233.7010.0!


Fitting of 3 drawers max. in addition to the standard drawer for table models 233.6110.0/233.7110.0!

**** Please indicate powder coating of frame and handles upon order placement (grey white RAL 9002, white aluminium RAL 9006 or grey aluminium RAL 9007).**

**Safe working load of medi-nets:
Totally: up to 50 kg
Drawer: up to 10 kg
Top/shelf: up to 20 kg**

Model No.

233.6010.0




Frame chromium-plated, 1 drawer, height 750 mm

233.7010.0



Frame powder coated, 1 drawer, height 750 mm

233.6110.0



Frame chromium-plated, 1 drawer, height 850 mm

233.7110.0



Frame powder coated, 1 drawer, height 850 mm

201.3000.0

1 additional drawer, with handle

Surcharges:

200.6600.0

for coloured frame acc. to colour chart (not possible for 233.6060.0)

200.6650.0

for coloured handle acc. to colour chart (not possible for 233.6060.0)

varimed® Accessories for medi-net general purpose and equipment tables

medi-net general purpose table 233.6010.0
with accessories for injecting location:
201.1170.0, 201.1610.0, 201.1700.0, 201.2600.0



medi-net general purpose table 233.6010.0
with accessories for ergometry
measuring workstations 201.1300.0
(w/o apparatus)



medi-net general purpose table 233.7010.0 +
201.0050.0 with accessories for home dialysis
201.1800.0, 201.1250.0 (w/o apparatus)



Cannula and syringe dispenser 201.1610.0



Console 201.1800.0



Infusion pole
201.1250.0



Cannula and syringe dispenser 201.1650.0



Ampule opener "Cupfix"
201.2600.0



"Medi-Müll" disposal set 201.2700.0




Refuse collector 201.2800.0



1 round stainless steel tray 201.1170.0



Drawer insert 201.2110.0 plastic




Drawer insert 201.2100.0 plastic

Model No.	Description	Model No.	Description	Model No.	Description
201.3000.0	1 drawer with handle, running rails, for subsequent fitting, for 580 mm wide medi-nets (internal dimensions: 510 x 413 x 90 mm, W x D x H)	201.0050.0	2 castors with locking device 75 mm Ø	201.1300.0	Cable guide fixing clamp 201.1910.0, diam. 16 mm
201.1170.0	1 round stainless steel tray diam. 320 mm, with fixing clamp 201.1900.0, diam. 10 mm	201.1250.0	Infusion pole € 2-section, with 4 hooks and fixing clamp 201.1920.0, diam. 18 mm, safe working load: 2 kg (per hook)	201.1610.0	Cannula and syringe dispenser grey, first in/first out system, with twin set incl. swab dispenser, requires fixing brackets 201.1700.0
201.1640.0	Cannula and syringe dispenser grey, first in/first out system, without right twin set, requires fixing brackets 201.1700.0	201.1650.0	Cannula and syringe dispenser grey, first in/last out system, with twin set incl. swab dispenser, requires fixing brackets 201.1700.0	201.1660.0	Cannula and syringe dispenser grey, first in/last out system, without twin set, requires fixing brackets 201.1700.0
201.1700.0	Fixing bracket 1 pair, for use with 201.1610.0, 201.1640.0, 201.1650.0 and 201.1660.0	201.1800.0	Console stainless steel with fixing brackets, for 580 mm wide medi-nets only	201.1900.0	Fixing clamp diam. 10 mm
201.2100.0	Drawer insert plastic, for 580 mm wide medi-nets only	201.2110.0	Drawer insert plastic, for 580 mm wide medi-nets only	201.1910.0	diam. 16 mm
201.2700.0	"Medi-Müll" disposal set 1,5 l 20 boxes made of polypropylene with red lid, with holder	201.2710.0	"Medi-Müll" spare boxes 1,5 l carton of 100 boxes, without holder	201.1920.0	diam. 18 mm
201.2800.0	Refuse collector with roll of bags, powder coated grey white	201.2850.0	Roll of spare bags for 201.2800.0	201.2600.0	Ampule opener "Cupfix" with fixing clamp
				201.2720.0	Holder as 201.2700.0, however without boxes. (Holders for different systems upon request)


Colours for powder coated finishes and upholsteries

Material description


Powder coated finishes



grey white RAL 9002



pure white RAL 9010



alabaster grey RAL 0007500*



light grey RAL 7035



anthracite grey RAL 7016




jet black RAL 9005



white aluminium RAL 9006



grey aluminium RAL 9007



colza yellow RAL 1009050*



signal yellow RAL 1003




pastel orange RAL 2003




light ivory RAL 1015




sahara beige RAL 0607030*



antique pink RAL 3014



heather violet RAL 4003




vermilion RAL 2002




light green RAL 6027




lemon green RAL 1008060*



leaf green RAL 1206050*



white blue RAL 2408015*




light blue RAL 2606030*



sapphire blue RAL 5003

Leatherette comfort



blossom white 84



gold yellow 86



mandarin 98



red carmine 92



raspberry red 103



caribbean blue 82



atlantic blue 83




topaz 105



opal green 100



apple green 94




lime 99

*RAL shades




milky coffee 102



natural brown 93



hazelnut 104



sand beige 101



stone grey 96




shale grey 97



titan grey 91

Leatherette classic




white 53



silver grey 55



ocean 52



black 57

Materials

Cabinet bodies, doors, claddings, and bases as well as are of sheet steel. Frames of round or rectangular steel tube.

Castors

All castors are resistant to abrasions.

Handles

of ergonomic design, distance of holes: 192 mm. The bowed bar handles are chemical resistant, easy to clean and to disinfect.

Powder coated finish

of high quality, multilayer gloss paint of extreme knock and scratch resistance. Easy to clean. Resistant to climatic conditions and disinfectants.

Standard finish of tubular steel frames: e.g. for couches: grey white RAL 9002, white aluminium RAL 9006, grey aluminium RAL 9007 or chromium-plated. Please refer to the different model descriptions.

Frames of couches, medi-net trolleys or stool bases are also available in different colours at extra charge. Please refer to our colour chart.

Standard finish of cabinet bodies is grey white RAL 9002; hinged doors grey white RAL 9002; standard finish of cabinet bases is jet black RAL 9005; handles are in grey white RAL 9002.

The SCHMITZ colours acc. to the colour chart are available without surcharge for handles.

Different finishes for cabinet bodies, hinges doors and bases acc. to our colour chart at extra charge.

Further special powder coatings acc. to RAL colour chart on request!

Upholstery

Foamed material on a support plate, with first-class, hard-wearing leatherette covers in two qualities.

Leatherette comfort:

18 colour tones

Wear-resisting material, with high resilience and flexibility providing excellent sitting and lying comfort. Grained surface. Available for an additional charge.

Please refer to our price list.

Leatherette classic:

4 colour tones

Hard-wearing, flexible material. Smooth surface.

Available at no extra charge.

We recommend:

Special cleaner for leatherette Model no. 2019878

500 ml spray bottle
(packaging unit 3 bottles)

Cleaning set for leatherette Model no. 2026973

1 x 500 ml spray bottle and a brush

Please refer to our separate operating manual for upholsteries for detailed information on cleaning and disinfection of upholsteries.

**SCHMITZ u. Söhne
GmbH & Co. KG**


P. O. Box 14 61
58734 Wickede (Ruhr)
Germany (postal address)
Zum Ostenfeld 29
58739 Wickede (Ruhr)
Germany (visiting address)
T +49 (0)2377 84 0
F +49 (0)2377 84 135
www.schmitz-soehne.com
export@schmitz-soehne.de

Technical service, hotline:
T +49 (0)2377 84 549
F +49 (0)2377 84 210
service@schmitz-soehne.de

All articles including correspond-
ing accessories marked with
the CE label meet the standards
for medical devices acc. to
MPG and the EC regulation
93/42/EEC.



The products may only be used
for their intended purpose of use.
Please also observe our operating
manual!




[f /schmitz.soehne](https://www.facebook.com/schmitz.soehne)



**Information on
our establishments
worldwide can be found
on our website.**

Gynaecology, Urology
and Proctology
The perfect equipment
to meet all requirements



vidan®
The Full HD Video Colposcope



Partura®
The SCHMITZ delivery table –
comfortable and safe delivery –
with focus on preventive
obstetrics

